# Propel Sinus Implant

## Patient Implant Information Leaflet

Always consult your physician regarding the information contained in this leaflet.

### 1. Device Description

#### a. Product Description

The PROPEL sinus implant releases a steroid (mometasone furoate) over time. The implant starts to break down and be absorbed once placed into the sinus. A delivery system is used to insert the implant.

#### b. Implant Material Composition

The Propel Sinus Implant is comprised of approximately 93 mg of polymers composed of lactide, glycolide, caprolactone, ethylene glycol and 370 mcg of mometasone furoate

### 2. Intended Use

The PROPEL<sup>®</sup> Sinus Implant is intended for use as a medical device in patients after sinus surgery. It is used to keep the ethmoid sinus open. It does this by spreading mucosal tissues, keeping the middle turbinate in place, stopping adhesions, and reducing swelling. The PROPEL<sup>®</sup> Sinus Implant reduces the need for some treatments after surgery, such as surgical adhesion lysis and/or use of oral steroids. Each PROPEL<sup>®</sup> Sinus Implant contains 370µg of mometasone furoate. This is released over time to reduce swelling.

### 3. Indications

The PROPEL<sup>®</sup> Sinus Implant is intended for use in patients age 18 years and older after sinus surgery. The implant keeps the sinus open, spreads apart nearby tissues, provides stability, prevents blockage, and reduces swelling. The Implant reduces the need for post-operative surgery and/or use of oral steroids.

### 4. Contraindications

The PROPEL sinus implant should not be used by the following patients:

- Patients that may be allergic to mometasone furoate or other steroids
- Patients with a known allergy to plastics, specifically lactide, glycolide or caprolactone copolymers.

### 5. Warnings

- The PROPEL Sinus Implant and delivery system are intended to be used once.

### 6. Precautions

- An allergic reaction may occur when the implant is placed into the sinus.
- In rare instances, sinus surgery, with and without sinus implants or packing, may present a risk of toxic shock syndrome (a rapid drop in blood pressure that can be fatal).
- This product has not been studied in Pregnant and Nursing Females.

### 7. Drug Interactions

No drug-drug reaction studies have been done with the implant. The risk of a drug-drug reaction with the implant is expected to be low.

### 8. Carcinogenicity, Genotoxicity, and Reproductive Toxicity

No long-term studies in animals have been done to assess the potential of the implant to cause cancer.

### 9. Lactation

There is currently no data that the steroid may affect and/or may be found in human milk.

### 10. Dosage and Administration

Each PROPEL implant contains 370µg of steroid which is released over time. The Device Lifetime for the Propel Sinus Implant is defined as 30 days. Talk to your doctor about any follow up that may be needed.

### 11. Potential Adverse Effects:

Potential side effects of the Propel sinus implant are thought to be like those found with other sinus stents, gels, or packing. Below is a list of the potential problems that may occur with use of the device:

- Early movement of implant or small implant pieces coming out of the nose
- Potential to swallow implant or implant pieces
- Crusting of implant causing, or worsening pain/pressure/headache
- Allergic reaction


Potential risks or side effects associated with intranasal mometasone furoate include:

- Nasal irritation
- Allergic reaction
- Nasal bleeding
- Infection (bacterial, fungal or viral) in the nose or throat
- Nasal burning
- Nasal dryness
- Increased risk of secondary infections due to bacteria, fungi or viruses
- Glaucoma/increase in eye pressure
- Cataracts/change in vision
- Headache
- Sore Throat

Potential risks or general side effects associated with steroids:

- Growth suppression
- Weakened immune system
- Allergic reactions
- Headache
- Nose bleeds
- Coughing
- Vomiting
- Fungal infection
- Glaucoma/increase in eye pressure
- Cataracts/changes in vision
- Joint pain
- Muscle pain

### 12. Instructions for completion of the Patient Implant Card

- Name of the Patient or Patient ID. To be filled by the healthcare institution/provider
- Date of implantation. To be filled by the healthcare institution or provider.
- Name and address of the healthcare institution/provider. To be filled by the healthcare institution/provider.

### 13. Symbols Used on the Patient Implant Card

	Name and Address of the implanting healthcare institution/provider
	Patient Name
	Date of procedure
	Medical Device
	Unique Device Identifier
	Manufacturer
	Lot number

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Intersect ENT Inc.