# PROPEL<sup>®</sup> Contour

(mometasone furoate implant, 370 µg)

## Instructions For Use

### CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE


**STERILE:** Sterilized by irradiation. Do not use if the package is open or damaged.  
**STORAGE:** The product should be stored at room temperature (approximately 25° C) with excursions permitted to 15-30° C.  
**SINGLE USE:** Product is supplied sterile and for single use only.

### PRODUCT DESCRIPTION

The PROPEL<sup>®</sup> Contour sinus implant provides sustained release of mometasone furoate via a bioabsorbable sinus implant. A delivery system is provided to insert the implant.


### Drug Component Description

The PROPEL Contour sinus implant contains mometasone furoate (active ingredient), a synthetic corticosteroid with anti-inflammatory activity. Mometasone furoate (MF) is a white to off-white powder. The chemical name is 9 $\alpha$ ,21-dichloro-11 $\beta$ ,17 $\alpha$ -dihydroxy-16 $\alpha$ -methylpregna-1,4-diene-3,20-dione 17-(2-furoate), with the empirical formula C<sub>27</sub>H<sub>30</sub>Cl<sub>2</sub>O<sub>6</sub>, and a molecular weight of 521.43 g/mol. Mometasone furoate is a hydrophobic drug that is practically insoluble in water. Mometasone furoate is stable under aqueous, acidic and oxidative conditions. Mometasone furoate can degrade under extreme basic, thermal and photolytic conditions. The chemical structure is shown below. The drug is embedded in a bioabsorbable polymer matrix containing poly-(DL-lactide-co-glycolide) and polyethylene glycol (inactive ingredients) which provides for gradual release of the drug.


Chemical structure of mometasone furoate

The inactive ingredients on the sinus implant are poly-(DL-lactide-co-glycolide) and polyethylene glycol. Poly-(DL-lactide-co-glycolide) is an amorphous biodegradable polymer. The chemical structure is shown below.

Chemical structure of poly-(DL-lactide-co-glycolide)

Polyethylene glycol is a hydrophilic polyether compound that is highly flexible. It is non-toxic and non-immunogenic. The chemical structure is shown below.

Chemical structure of polyethylene glycol

### Implant Component Description

The PROPEL Contour sinus implant is comprised of a bioabsorbable co-polymer, poly(L-lactide-co-glycolide), PLG and poly(L-lactide-co- $\epsilon$ -caprolactone).

The implant is bioabsorbable and is designed to accommodate the size and variability of the surgically enlarged frontal sinus ostium. Once inserted, the implant is designed to be self-retaining against the mucosa of the surgically enlarged sinus ostium in order to maintain sinus patency and deliver drug to the mucosa. The implant should be inserted by a physician under endoscopic visualization. A delivery system is provided to access the frontal sinus ostium and deploy the implant. A loading tool and funnel are optional accessories provided to assist in loading of the implant into the delivery system.

### INDICATIONS AND INTENDED USE

The PROPEL Contour sinus implant is intended for use in patients  $\geq$  18 years of age with chronic rhinosinusitis following sinus surgery to maintain patency of the frontal sinus ostia and to locally deliver mometasone furoate to the sinus mucosa. The PROPEL Contour sinus implant separates/dilates mucosal tissues, prevents obstruction by adhesions/scarring, and reduces edema. The implant reduces the need for post-operative intervention such as surgical adhesion lysis and/or use of oral steroids.

### CONTRAINDICATIONS:

- The use of the PROPEL Contour sinus implant is contraindicated in the following patients:
- Patients with suspected or confirmed hypersensitivity and/or intolerance to mometasone furoate.
  - Patients with a known hypersensitivity and/or intolerance to lactide, glycolide or caprolactone copolymers.

### WARNINGS

- The PROPEL Contour sinus implant is designed for single patient use only. Do not reuse, reprocess, or re-sterilize. Reuse, reprocessing, and re-sterilization may compromise the structural integrity of the device and/or lead to device failure that may result in patient injury. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury or illness of the patient.
- Do not use if the package is open or damaged.

### PRECAUTIONS

- Special care should be taken to avoid bending, twisting or damaging the implant.
- The implant is not designed to be modified by the physician.
- The implant is not intended to be compressed and loaded into the delivery system more than two times.
- The implant must be placed under endoscopic visualization.
- The implant exhibits no antimicrobial properties.
- Foreign body reaction may occur as is possible with most surgical adjuncts.
- In rare instances, the physiochemical condition associated with sinus surgery, both with and without sinus implants or packing, may present a risk of toxic shock syndrome (TSS).
- Pediatric Use: The safety and effectiveness of the implant in pediatric patients have not been established.
- Pregnancy and Nursing Females: The safety and effectiveness of the implant in pregnant or nursing females have not been established.
- The safety and performance of the implant in immunodeficient patients, patients with active systemic bacterial, fungal, viral infection, and patients with insulin dependent diabetes have not been established.

### DRUG INFORMATION

**MECHANISM OF ACTION:** Corticosteroids have been shown to have a wide range of effects on multiple cell types (e.g., mast cells, eosinophils, neutrophils, macrophages, and lymphocytes) and mediators (e.g., histamine, eicosanoids, leukotrienes, and cytokines) involved in inflammation. The precise mechanism behind the anti-inflammatory properties of the eluted mometasone furoate is not known.

**PHARMACOKINETICS:** The PROPEL Sinus Implant underwent pharmacokinetic testing. Following bilateral drug-eluting PROPEL implant placement after endoscopic sinus surgery for chronic rhinosinusitis and subsequent weekly morning blood sampling for 4 weeks in 5 adult patients, plasma mometasone furoate concentrations were not quantifiable at any time point. Mean cortisol concentrations were within normal limits.

### DRUG INTERACTIONS

No drug-drug interaction studies have been conducted with the implant.

### CARCINOGENICITY, GENOTOXICITY AND REPRODUCTIVE TOXICITY

No long term studies in animals have been performed to evaluate the carcinogenic potential of the implant.

### PREGNANCY

There have been no controlled studies in pregnant women using the PROPEL Contour sinus implant. The PROPEL Contour sinus implant should be used during pregnancy only if the potential benefits justify the potential risk.

### LACTATION


It is not known if mometasone furoate is excreted in human milk. Because other corticosteroids are excreted in human milk, the implant should be used only if the potential benefits justify the potential risk.

### DOSAGE AND ADMINISTRATION

Each PROPEL Contour sinus implant contains 370 µg of mometasone furoate which is gradually released over time.


### DIRECTIONS FOR USE

1. Remove the implant and delivery system from its protective packaging using sterile technique. Inspect for any obvious damage.
2. Prepare the delivery system.
  - a. Depress and hold the plunger.
  - b. Using your fingers, gently bend the tip of the applicator to the desired angle while supporting the tubing with your fingers. The useable range of the delivery system is from 70 to 110 degrees. Note: It may be necessary to bend the tip past the desired angle to achieve the appropriate angle.
  - c. Release the plunger.

Step 2a


Step 2b

Step 2b

3. The implant must be compressed and loaded into the tip of the delivery system prior to use. Note: The compressed implant may be loaded directly into the tip of the delivery system or with the use of the optional funnel and loading tool accessories.

- a. For implant loading with the use of the funnel, place the funnel onto the tip of the delivery system.
- b. Carefully remove the top of the implant holder, being careful not to grasp or dislodge the implant.
- c. Gently lift the implant off of the base of the implant holder.
- d. Grasp the implant between the fingers of both hands and gently compress the implant.
- e. Insert compressed implant into either the attached funnel (if applicable) or directly into to the distal tip of the delivery system.
- f. If the funnel is attached, use the loading tool to push the implant past the opening of the funnel. If the funnel is not attached, continue to push the compressed implant directly into the delivery system until it is flush with the distal tip.
- g. If the funnel is attached, carefully remove the funnel, taking care not to dislodge the implant from the tip of the delivery system. If the implant begins to withdraw from the tip during funnel removal, replace the funnel and gently squeeze the tip of the delivery system to hold implant in place.
- h. The implant may be compressed and loaded into the delivery system tip up to two times. **CAUTION:** Do not leave the implant in the crimped state for more than three minutes prior to placement.

Step 3b

Step 3d

Step 3f

Step 3g

Step 5b

4. For adequate visualization, ensure hemostasis in operated sinus ostia and/or cavities prior to insertion. Advance the delivery system into the sinus cavity using endoscopic visualization.
5. To insert the implant in the sinus ostium:
  - a. Ensure the delivery system is oriented so the distal tip is curved appropriately for the targeted patient anatomy and advance the distal tip of the delivery system into the sinus ostium.
  - b. Note: The markings (shown above) on the distal end of the delivery system may be used to aid in placement (for reference only).
  - c. Deploy the implant halfway by depressing the plunger until the distal edge of the coiled pusher aligns with the green to clear transition. At this point, the implant is approximately halfway deployed.

# PROPEL® Contour

(mometasone furoate implant, 370 µg)

## Instructions For Use

- d. Slightly withdraw the delivery system to visualize the waist of the implant (which is denoted by the strut crossing-points).
  - e. Align the waist of the implant with the ostium.
  - f. While keeping the waist of the implant aligned with the ostium, continue withdrawing the delivery system from the sinus ostium while simultaneously continuing to depress the plunger.
  - g. Visualize the waist of the implant and ensure it is aligned with the ostium before completely depressing the plunger and fully deploying the implant.
6. Confirm final placement by endoscopic visualization. Confirm the waist of the implant is aligned with the sinus ostium. Confirm the implant is well apposed to the tissue to maximize drug delivery. To adjust the position of the implant, use standard surgical instruments.

### POST-OPERATIVE CARE

- As part of routine post-operative care, frequent use of saline sprays, rinses or irrigations is recommended to keep the implant moist.
- Routine debridement may be performed as part of the usual post-operative care.
- The implant may be removed at the discretion of the physician by use of suction, forceps or other surgical instruments.

### CLINICAL TRIALS

The safety and efficacy of the PROPEL Contour sinus implant when placed in the frontal sinus ostia (FSO) following endoscopic sinus surgery (ESS) with traditional instrumentation, balloon dilation, or a combination of both in patients with chronic rhinosinusitis were assessed in 2 prospective clinical trials conducted in the United States.

The EXCEED study was a prospective, single-arm, open-label, feasibility trial designed to assess the performance, safety and initial signals of efficacy of the PROPEL Contour sinus implant when used in patients with chronic rhinosinusitis following ESS in the office and operating room settings. Fifteen of 16 patients enrolled at 2 study sites were eligible to participate. The EXCEED study results showed an overall implant delivery success rate of 97.8% (100% in frontal sinuses; 95.2% in maxillary sinuses), demonstrating that the primary performance endpoint was met. No serious unanticipated AEs were reported; therefore, the study met the primary safety endpoint. In addition, no AEs were related to placement of the PROPEL Contour sinus implant. Observed endoscopic improvements at Day 90 included a 46.8% increase in ostial patency, a 75.0% reduction in polypoid edema, a 19% reduction in grade 2 and 3 adhesion/scarring, and a 38.1 mm decrease in inflammation score. The mean SNOT-22 score at Day 90 showed a treatment effect size of 1.09, representing a large health benefit. None of the patients required oral steroids for sinus obstruction, and only 2 received surgical interventions (i.e., polypectomy; removal of left frontal scar) through the duration of the study.

The PROGRESS study was a prospective, randomized, blinded, controlled study that enrolled 80 patients across 12 study centers. The study utilized an intra-patient control design to assess the safety and efficacy of the PROPEL Contour sinus implant when placed following surgery on 1 sinus side compared to surgery alone on the contralateral side. The primary efficacy endpoint was the reduction in need for post-operative interventions at Day 30, as determined by an independent, blinded, sinus surgeon based on video-endoscopy review. Post-operative intervention was a composite endpoint that included surgical intervention required to debride obstructive adhesions or scar tissue formation (grade 2 or 3) in the frontal sinus ostia (FSO), and/or oral steroid intervention warranted to resolve recurrent inflammation or polypoid edema in the frontal recess/FSO. Secondary efficacy endpoints of frequency and severity of adhesion/scarring, polypoid edema and inflammation were determined endoscopically by an independent reviewer and clinical investigators at the study centers. The safety measures of adverse events and serious adverse events were recorded throughout the 90-day follow-up period.

The PROGRESS study demonstrated that PROPEL Contour sinus implant placement in the FSO was successful (100% implant delivery), safe, and effective in significantly reducing the need for post-operative interventions in the FSO at Day 30 (p=0.0023), as judged by an independent blinded reviewer when no other post-operative corticosteroids were administered for 14 days (see table below). Holm's step-down procedure used to control for family-wise type I error rate on 6 prespecified secondary efficacy endpoints showed statistically significant reductions in the need for postoperative interventions (p=0.0039), need for surgical interventions (p=0.0156), degree of inflammation (p=0.0005), and rate of occlusion/restenosis (p<0.0001) as observed by clinical investigators at Day 30. Through Day 90, the rate of occlusion/restenosis was significantly (p=0.0055) lower on the treatment compared to the control sides, representing a relative reduction of 42.9%.

	Treatment (Tx) (n=80)	Control (Ctrl) (n=80)	p-value
<b>PRIMARY EFFICACY RESULTS<sup>a,†</sup></b>			
Need for Post-Operative Interventions, N (%)	7 (11.5%) (n=61) <sup>§</sup>	20 (32.8%) (n=61) <sup>§</sup>	0.0023
<b>SECONDARY EFFICACY RESULTS<sup>b,†</sup></b>			
Need for Post-Operative interventions, N (%)	12 (16.0%) (n=75) <sup>§</sup>	25 (33.3%) (n=75) <sup>§</sup>	0.0039
Need for Surgical Interventions, N (%)	3 (4.0%) (n=75) <sup>§</sup>	11 (14.7%) (n=75) <sup>§</sup>	0.0156
Inflammation (100-VAS, mm), Mean (SD)	23.1 (24.23) (n=79) <sup>§</sup>	35.6 (31.12) (n=77) <sup>§</sup>	0.0005
Occlusion/Restenosis, N (%)	10 (13.3%) (n=75) <sup>§</sup>	27 (36.0%) (n=75) <sup>§</sup>	<0.0001

<sup>a</sup>All eighty patients returned for the Day 30 visit and had their endoscopy recorded for grading by independent reviewer; however, data were considered missing if the independent reviewer could not grade a video due to sub-optimal video quality or inadequate imaging of the relevant anatomy. Inadequate imaging of the relevant anatomy can occur when presence of significant edema or an adhesion prevents access into the frontal sinus. Since the planned statistical test (McNemar's test of correlated proportions) requires patients with an observed pair of outcomes, 19 patients could not be included in the test. McNemar's exact binomial test was employed to obtain the 2-sided p-value at alpha level of 0.05 for the primary efficacy endpoint and other categorical efficacy endpoints; T-tests were performed for all continuous efficacy data on the side-to-side difference in scores

<sup>†</sup>p-values adjusted using a Holm's step-down method to control for familywise type I error

<sup>b</sup>Determined at Day 30 by the independent reviewer based on video-endoscopy review

<sup>†</sup>Determined at Day 30 by clinical investigators

<sup>§</sup>Number of patients with evaluable sinuses

SD=Standard Deviation, VAS=Visual Analog Scale

Please refer to <https://clinicaltrials.gov/> for additional details on the EXCEED (NCT02228720) and PROGRESS (NCT02266810) clinical studies.

### ADVERSE EVENTS

#### OBSERVED ADVERSE EVENTS

PROPEL Contour sinus implant is a drug-releasing sinus implant that is smaller in size than the PROPEL Mini, to accommodate the size and variability of the enlarged sinus ostium. Two prospective clinical trials (EXCEED and PROGRESS) conducted in the United States studied 95 patients and 125 PROPEL Contour sinus implants. In the EXCEED study, no implant-related serious adverse events were reported, resulting in a 0% incidence rate of implant-related adverse events. No patients withdrew due to an adverse event, and no deaths occurred.

In the Nova cohort (with PROPEL Contour) of the PROGRESS study with 80 patients, no implant-related serious adverse events or adverse events were observed, resulting in a 0% overall incidence rate of implant-related adverse events. Three adverse events (headache, epistaxis, acute sinusitis) were judged by clinical investigators as having an indeterminate relationship to the implant. All three events resolved without sequelae. No patients withdrew due to an adverse event and no deaths occurred. Adverse events (regardless of relationship to implant), which were reported in ≥ 2% of patients in the PROGRESS study are summarized in the table below.

Adverse Events (n=80) – PROGRESS study (Nova Cohort)	
Adverse Event Type	Patients Reporting (%)
Acute Sinusitis	20.0
Asthma	7.5
Headache	6.3
Chronic Sinusitis	5.0
Upper Respiratory Tract Infection	5.0
Fungal Sinusitis	2.5

Adverse Events (n=80) – PROGRESS study (Nova Cohort)	
Adverse Event Type	Patients Reporting (%)
Nasopharyngitis	2.5
Nausea	2.5
Neck Pain	2.5
Sinus Headache	2.5
Streptococcal Pharyngitis	2.5

#### POTENTIAL ADVERSE EVENTS

Risks associated with the use of the PROPEL Contour sinus implant are anticipated to be similar to those experienced by patients who undergo endoscopic sinus surgery, and placement of sinus implants or packing.

Additional risks potentially associated with placement of the PROPEL Contour sinus implant may include, but not be limited to the following.

- Premature displacement of implant or implant fragments
- Swallowing of implant or implant fragments
- Pain/pressure/headache due to the adherence of crusting to, or presence of the implant
- Aspiration of small implant fragments
- Foreign body response, including formation of granulation tissue
- Cerebrospinal Fluid (CSF) leak from surgical procedure.

Risks potentially associated with intranasal mometasone furoate may include, but not be limited to the following.

- Nasal irritation
- Hypersensitivity reaction
- Intranasal bleeding (epistaxis)
- Localized infection (bacterial, fungal or viral) in the nose or pharynx
- Nasal burning
- Nasal dryness
- Susceptibility to secondary infections due to bacteria, fungi or viruses
- Glaucoma/elevation of intraocular pressure
- Cataracts/change in lens opacities
- Headache
- Pharyngitis

Risks potentially associated with steroids may include, but not be limited to the following.

- Alteration of the HPA axis including growth suppression
- Immunosuppression
- Hypersensitivity reactions
- Headache
- Epistaxis
- Coughing
- Vomiting
- Candidiasis
- Glaucoma/elevation in intraocular pressure
- Cataracts/changes in lens opacities
- Arthralgia
- Myalgia

There may be other potential adverse effects that may occur which are currently unforeseen.

Symbols Used on Product Labeling									
	Catalogue Number		Do Not Re-Use		Use By		Room Temperature		Customer Service Number
	Batch Code		Sterilized using Irradiation		Caution		Do Not Use If Package Is Damaged		
	Do Not Resterilize		Consult Instruction For Use		Authorized Representative In The European Community		Manufacturer		

#### Product Information Disclosure

Intersect ENT, Inc. has exercised reasonable care in the manufacture of this product. Intersect ENT, Inc. excludes all warranties, whether expressed or implied, by operation of law or otherwise, including but not limited to, any implied warranties of merchantability or fitness, since handling and storage of this product, as well as factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Intersect ENT, Inc.'s control, directly affect this product and the results obtained from its use. Intersect ENT, Inc. shall not be liable for any incidental or consequential loss, damage or expense, directly or indirectly arising from the use of this product. Intersect ENT, Inc. neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this product.

The above-identified implants, delivery systems, and/or the use of the above-identified implants/delivery systems in a method may be covered by one or more U.S. and/or foreign patents, found at <https://intersectent.com/patents>.

© 2021 Intersect ENT Inc. All rights reserved. INTERSECT ENT® and PROPEL® are registered trademarks of Intersect ENT, Inc.

Intersect ENT Inc.  
1555 Adams Drive  
Menlo Park, CA 94025 USA  
[www.intersectENT.com](http://www.intersectENT.com)

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© 2021 Intersect ENT Inc. Med enerett. INTERSECT ENT® og PROPEL® er registrerte varemerker for Intersect ENT, Inc.